






CIRCULAR STAPLER



For different calibers the corresponding staple legs will be triggered into the right figure of »B« for the different tissue thicknesses.



Article	Order No.	Colour of Protector	Number of Staples	Outer Diameter*	Inner Diameter*	Leg Height*	Sales Unit
WH-Y-17	WH-Y-17	white	16	17.0	10.0	3,8	3 pcs
WH-Y-21	WH-Y-21	yellow	20	21.0	14.0	4,5	3 pcs
WH-Y-24	WH-Y-24	grey	18	24.0	15.0	4,5	3 pcs
WH-Y-26	WH-Y-26	red	20	26.0	17.0	4,5	3 pcs
WH-Y-29	WH-Y-29	blue	24	29.0	20.0	4,8	3 pcs
WH-Y-32	WH-Y-32	green	26	32.0	22.0	5,0	3 pcs
WH-Y-34	WH-Y-34	purple	28	34.0	24.0	5,0	3 pcs