



Skin Stapler



Remover

Einweg-Skin Stapler zum Verschließen von Wunden und chirurgischen Inzisionen. Mit dem Klammerentferner lassen sich die Hautklammern einfach entfernen.

Material: Die Klammern bestehen aus rostfreiem Stahl 00Cr18Ni14Mo3 gemäß dem GB4234-2003. Der Griff besteht aus ABS gemäß GB12671. Der Klammerentferner besteht aus ABS und rostfreiem Stahl 1Cr18Ni9.

Leistungsmerkmale: Steril, gute Bioverträglichkeit, gute Metallstabilität, schnelle Durchführung der chirurgischen Naht, klare Sicht auf die Wundränder, einsetzbar bei allen Schnittführungen.

Bestell-Nr.: Skin Stapler: PH-W-35-12

Remover: PC-01-12

pro VE: 12 Stk.



SKIN STAPLER & REMOVER



Skin Stapler



Remover

Single-use Skin Stapler for wound closure and surgical incisions. Using the Skin Stapler Remover, the skin staples can be easily removed.

Material: The staples consist of stainless steel 00Cr18Ni-14Mo3 according to the GB4234-2003. The handhold is made of ABS according to GB12671. The Skin Stapler Remover consists of ABS and stainless steel 1Cr18Ni9.

Performance features: sterile, excellent biocompatibility, very good metal stability, quick realization of the surgical suture, clear view to the wound edges, suitable for all cross sections.

Reference-No.: Skin Stapler: PH-W-35-12 / Packing unit 12 pcs

Remover: PC-01-12 / Packing unit 12 pcs

CE0044