


Disposable Linear Cutter Stapler for Endoscopic Use



- A) Pin
- B) Unload/Unlock Button
- C) Shaft
- D) Articulating lever
- E) Rotation collar
- F) Green button
- G) Black return knob
- H) Handle
- I) Single use loading unit



- J) Load alignment indicators
- K) Shipping wedge
- L) Anvil
- M) Tissue stop
- N) End of staple line
- O) End of cut line
- P) Increment markings
- Q) Lower clamp button
- R) Staple cartridge

Disposable Linear Cutter Stapler for Endoscopic Use

The Disposable Flexible Linear Cutter Stapler has applications in abdominal, gynecologic, pediatric and thoracic surgery for resection, transection and creation of anastomosis.

Article	Instrument Order No.	Cartridge Order No. (for both stapler sizes)	Leg Height	Number of Staples	Cartridge color
Endocutter Stapler Standard Model	ZNQ-A	ZN-Z-45-2,5	2,5 mm	66	white
		ZN-Z-45-3,5	3,5 mm	66	blue
		ZN-Z-45-4,8	4,8 mm	66	green
Endocutter Stapler XL Model	ZNQ-C	ZN-Z-60-2.5	2,5 mm	90	white
		ZN-Z-60-3.5	3,5 mm	90	blue
		ZN-Z-60-4.8	4,8 mm	90	green

CE 0197