


LINEAR CUTTER



Disposable Stapler with Titanium Staples

1. Cartridge
2. Cartridge Housing
3. Anvil
4. Protector
5. Indicator
6. Handle
7. Refire

Available sizes: 55 mm / 75 mm / 100 mm

Article	Instrument Order No.	Cartridge Order No.	Tissue Thickness*	Leg Height*	Number of Staples	Cartridge color	Sales Unit
QHB-55	QHB55	QCB55	1.5	3.85	56	blue	3 pcs
QHH-55	QHH55	QCH55	2.0	4.50	56	green	3 pcs
QHB-75	QHB75	QCB75	1.5	3.85	76	blue	3 pcs
QHH-75	QHH75	QCH75	2.0	4.50	76	green	3 pcs
QHB-100	QHB100	QCB100	1.5	3.85	100	blue	3 pcs
QHH-100	QHH100	QCH100	2.0	4.50	100	green	3 pcs

*mm

