

SKIN STAPLER & REMOVER



Skin Stapler
CE0120



Remover
CE0197

Single-use Skin Stapler for wound closure and surgical incisions.
Using the Skin Stapler Remover the skin staples can be easily removed.

Material:

The staples consist of stainless steel according to AISI 316L.
The handhold is made of ABS according to CAS 9003-56-9. The Skin Stapler Remover consists of ABS and stainless steel 1Cr18Ni9.

Performance features:

sterile, excellent biocompatibility, very good metal stability, quick realization of the surgical suture, clear view to the wound edges, suitable for all cross sections.

Order No.	Size of staples	Sales Unit
Skin Stapler: 35W	6.4 × 4.0 mm	50 pcs
Remover: SRK-T-2		50 pcs