

SURGISPON – Absorbable Haemostatic Gelatin Sponge



SURGISPON Tampon anal/vaginal

SURGISPON is a quick and effective hemostatic which stops bleeding fast, generally within 2 minutes. It reduces intraoperative bleeding and absorbs 40 – 50 times of its own weight of water/oxalate blood.

SURGISPON is non-toxic, non-illergenic, non-immunogenic and non-pyrogenic. SURGISPON is already gamma-sterilized and available in a convenient, ready to tear sterile blister/envelope. Unlike many haemostatic agents, it does not require special storage conditions.

To achieve fast hemostasis in:

- Haemorrhoidectomies, anal fistula and fissures
- Curettage and Hysterectomies

SURGISPON reduces significantly postoperative pain in 83% of the cases. It may be used in combination with antibiotics, Thrombin, chemo therapeutics without reduction of the haemostatic effect.

SURGISPON may provide vast economical benefits in recto-vaginal surgery. The gelatin sponge application can save up to 4 days of hospital stay for each rectal patient and valuable time for surgeons. Surgispon is easily expelled at first bowel movement, requiring no additional nursing intervention. If left in situ, it will be absorbed.

Order No.	Product	Dimension (mm)	Sales Unit
SSP-8030	Tampon	80 × 30 Ø	24 pcs

CE 0434