



Skin Stapler



Remover

Einweg-Skin Stapler zum Verschließen von Wunden und chirurgischen Inzisionen. Mit dem Klammerentferner lassen sich die Hautklammern einfach entfernen.

Material:

Die Klammern bestehen aus rostfreiem Stahl gemäß AISI 316L.
Der Griff besteht aus ABS gemäß CAS 9003-56-9. Der Klammerentferner besteht aus ABS und rostfreiem Stahl 1Cr18Ni9.

Leistungsmerkmale:

Steril, gute Bioverträglichkeit, gute Metallstabilität, schnelle Durchführung der chirurgischen Naht, klare Sicht auf die Wundränder, einsetzbar bei allen Schnittführungen.

Bestell-Nr.:	Klammergröße	pro VKE:
Skin Stapler: 35W	6,4 × 4,0 mm	10 Stk.
Remover: SRK-T-2-10		10 Stk.



SKIN STAPLER & REMOVER



Skin Stapler



Remover

**Single-use Skin Stapler for wound closure and surgical incisions.
Using the Skin Stapler Remover the skin staples can be easily removed.**

Material:

The staples consist of stainless steel according to AISI 316L.
The handhold is made of ABS according to CAS 9003-56-9. The Skin Stapler Remover consists of ABS and stainless steel 1Cr18Ni9.

Performance features:

sterile, excellent biocompatibility, very good metal stability, quick realization of the surgical suture, clear view to the wound edges, suitable for all cross sections.

Product Code:	Size of staples	Pieces per unit
Skin Stapler: 35W	6.4 × 4.0 mm	10
Remover: SRK-T-2-10		10

